



Quick Thermal Conductivity Meter

QTM-500

**Probe Sensors to Any of Your Application Needs !
Thermal Conductivity of All Kinds of Sample Materials
by Just One TC Meter!**



KYOTO ELECTRONICS

Quick Thermal Conductivity Meter

All kinds of sample can be measured by two different types of sensor up to your needs!

QTM-500



QTM-500 does quick and easy measurement of all kinds and types of sample materials. Stick sensor probe on sample surface of thermal equilibrium, and press START key. The measurement results will appear on display in 60 seconds.


Probe

Principle of measurement


The probe consists of single heater wire and thermocouple. When constant electric power (energy) is given to the heater, the temperature of the wire will rise in exponential progression. Temperature rising curve is plotted in linear line in below figure with time axis scaled in logarithm. The angle of this line increases if the sample has less thermal conductivity, and decreases if it has higher TC. Therefore, TC of a sample can be determined from the angle of the rising temperature graphic line.

$$\lambda = \frac{q \cdot \ln(t_2/t_1)}{4 \pi (T_2 - T_1)}$$

λ ; thermal conductivity of sample [W/mK]
 q ; generated heat per unit length of sample/time [W/m]
 t_1, t_2 ; measured time length [sec]
 T_1, T_2 ; Temperature at t_1, t_2 [K]



● Box type probe (standard) PD-11



● Insulation damp - proof probe (option)

For measurement of hydrous or electrically conductive material like foods or crude concrete, etc.


PD-13



Examples of LC Display

● Start measurement

When START button is pressed, the heater is supplied with constant power by heater current. The graphic curve shows realtime thermal progression of sample surface (heater temperature).




● Measurement result

The measurement results will be shown on display immediately after the measurement is finished.



● Measurement results (Temperature graphic for time log)

Temperature change can be confirmed by linear graphic line for time in logarithm.



Other Options

Printer DP-500

Sample No., TC and temperature are printed out.



Powder measuring case QTM-PA1

Powder sample can be measured by combination of this case and PD-11 probe.

SOFT-QTM5EW

Software for Measuring Thermal Conductivity of Thin Sheet (Optional parts)

Typical example of measurement

The QTM-500 Thermal Conductivity Meter and a personal computer with Windows® 95 were connected, and SOFT-QTM5EW software was installed. Then, measurements for thermal conductivity were performed on the following samples; Homogeneous material in film, sheet or thin board form, of those sample like rubber, plastics, ceramics, paper, textile or wood.

Thickness of sample

30 μ m to 10mm thick material in sheet form

Measuring range

0.035 ~ 5.0w/m K


Measuring equipment

CD-ROM _____ 1pce
 Operating manual _____ 1copy
 PC connecting cable for DOS-V _____ 1pce


● Main Menu for Measurement



● Measurement of Reference Material



● Measurement of Thin Sheet



■Specification QTM-500

Type and model name	QTM-500 Quick Thermal Conductivity Meter
Measuring method	Hot Wire method
Measuring range	0.023 to 12W/mk
Precision	±5% reading value per reference plate
Reproducibility	±3% reading value per reference plate
Temperature	-10 to 200°C (Thermal bath is necessary for measurement at different temperature than room)
Measuring time	Standard 60sec (specimen must be in equilibrium)
Sensor	PD-11 Box Probe Constantan heater wire and chromel-alumel thermocouple
Heater current precision	±0.05% of setup value
Display	30 characters 7 lines LCD with back light Display: Thermal conductivity: 0. XXXX~XXX. XW/mK Measuring temperature: -100 to 1000°C Guiding message for measurement
Minimum sample required	approx. 100W×50L×20 thick mm or more
External output	RS - 232C 2 channels Thermocouple electro-motive force 1 channel

Ambient condition	Temperature: 5 to 35°C Humidity: below 85%RH																							
Power source	AC100 to 240V	50/60Hz																						
Power consumption	60W																							
Dimension	Main unit: 300W×475DX175Hmm Box probe: 100W×50DX100Hmm																							
Weight	9kg																							
Supplied parts	<table border="0"> <tr> <td>1. PD-11 Box probe</td> <td>1pce</td> </tr> <tr> <td>2. Probe constant card</td> <td>1pce</td> </tr> <tr> <td>3. Power cord with ground wire</td> <td>1pce</td> </tr> <tr> <td>4. Power fuse</td> <td>2pce</td> </tr> <tr> <td>5. Reference plate:</td> <td></td> </tr> <tr> <td> R1-2 Clear quartz in box</td> <td>1pce</td> </tr> <tr> <td> R2-2 silicone rubber in box</td> <td>1pce</td> </tr> <tr> <td> R3-2 polyethylene form in box</td> <td>1pce</td> </tr> <tr> <td>6. Cooling plate (aluminum)</td> <td>1pce</td> </tr> <tr> <td>7. Brush</td> <td>1pce</td> </tr> <tr> <td>8. Operating manual</td> <td>1copy</td> </tr> </table>		1. PD-11 Box probe	1pce	2. Probe constant card	1pce	3. Power cord with ground wire	1pce	4. Power fuse	2pce	5. Reference plate:		R1-2 Clear quartz in box	1pce	R2-2 silicone rubber in box	1pce	R3-2 polyethylene form in box	1pce	6. Cooling plate (aluminum)	1pce	7. Brush	1pce	8. Operating manual	1copy
1. PD-11 Box probe	1pce																							
2. Probe constant card	1pce																							
3. Power cord with ground wire	1pce																							
4. Power fuse	2pce																							
5. Reference plate:																								
R1-2 Clear quartz in box	1pce																							
R2-2 silicone rubber in box	1pce																							
R3-2 polyethylene form in box	1pce																							
6. Cooling plate (aluminum)	1pce																							
7. Brush	1pce																							
8. Operating manual	1copy																							



KYOTO ELECTRONICS MANUFACTURING CO.,LTD.

Orerseas Division : 8-3 Niban-cho Chiyoda-ku TOKYO 102-0084,
JAPAN Fax : +81-3-3237-0537,Phone : +81-3-3239-7331

URL: <http://www.kyoto-kem.com>

Distributed by